



Product Main Number 1000021316



Shaped pad, ABENA San, 11/XXL, purple colourcode, Premium

- ✓ Longer pad to suit a range of body shapes
- ✓ High performance core with Top-Dry Technology
- ✓ Visible and invisible leakage barriers for extra protection
- ✓ Wetness indicator with graduation scale
- ✓ Support healthy skin balance



Product Description

ABENA San 7/XXL and 11/XXL are ABENA's range of anatomically shaped pads designed for bariatric use. The longer pad provides greater security for all body shapes. The high performance core with Top-Dry Technology ensures optimal skin care by keeping the exposed area comfortable and dry even after multiple voids. The fully breathable pads are fitted with standing leakage barriers combined with built-in cross barriers at the rear. When combined, these functionalities help protect the skin and maintain a natural pH skin balance. The Odor control system minimizes unpleasant smells. The wetness indicator with graduation scale guides the user and caretaker to determine the right absorption level and product changes. The dermatologically tested products have been given the Nordic Swan Ecolabel. These bladder protection pads should always be used with special fixation underwear, like our ABENA Fix Pants Super 4-6XL products, both to ensure a tight fit and to reduce the risk of leakage. ABENA's two-piece systems, consisting of a fixation pants and an incontinence pad, are a cost-effective solution, as they use less raw materials, thereby having less impact on natural resources and the environment.



Specifications

Base name	Shaped pad
Brand	ABENA
Sub-Brand	San
Size	11/XXL
Color	White
Features	Purple colourcode
Material	ECF fluff, PE, breathable textile backsheet (BTBS), elastane, nonwoven, PP, resin, SAP
Length/depth	84 cm
Width	46 cm
Tolerance	+/- 10%
Absorption Capacity	3400 ml
Gender	Unisex
Quality level	Premium
Certificates	CE. Nordic Swan Ecolabel. MD. UKCA.
CE Class (Medical Devices)	Class I
Product or test standards	ISO 11948-1, ISO 13485, ISO 14001, ISO 9001
Directives, regulations and acts	MDR (EU) 2017/745

Shelf Life	5 years
Storage Instructions	Store dry, room temperature and no direct sunlight.
Product Disposal Instructions	Can be disposed of with normal household waste sorted according to local regulations.
Packaging Disposal Instructions	Can be recycled or incinerated.

Instructions for use/application

Anatomically shaped pads for moderate to heavy incontinence specially designed for bariatric use thanks to the additional pad length. Easy product recognition by color codes, names, and numbering systems. These bariatric pads should always be used with special fixation underwear, like our ABE-NA fix Pants Super 4-6XL. That way the pad is kept in place, ensuring a tight fit and reducing the risk of leakages. Pull the fixation pants half-way up the legs. Fold the pad lengthways to create a bowl shape. Pass the pad between the legs from front to back. Position the upper rear end of the pad at the tailbone. Fan out the back of the pad horizontally to create a channel shape between the legs. Use the wetness indicator to check if the product is applied accurately in the center. Fan out the pad horizontally at the front. Pull the fixation pants up into position. See the application guide in our brochure for detailed information.

Packaging data

Unit	Contains	Length	Width	Height	EAN
cII	4 pck	57.6 cm	39.6 cm	36.7 cm	5713571000755
pck	16 pcs	38 cm	18 cm	26 cm	5713571000748



Nordic Swan Ecolabel is an official Scandinavian environmental label. Products with the Nordic Swan Ecolabel meet stringent environmental requirements with regard to raw material, production, use, documentation, and safety.



Regulation (EU) 2017/745 is a regulation of the European Union on the clinical investigation and sale of medical devices for human use. It repeals Directive 93/42/EEC, which concerns medical devices, and Directive 90/385/EEC, which concerns active implantable medical devices, on 26 May 2021.



The CE mark guarantees that a product is safe to use and complies with all safety precautions. CE stands for Conformité Européenne (European Conformity) and is mainly found on electronic equipment, safety equipment, construction products, and medical devices.



The UKCA mark guarantees that a product is safe to use and complies with all safety precautions in line with product standards and legislation of the UK.



Your guarantee that our products have been tested for ingredients which can potentially cause skin irritation. Tests conducted by an independent third-party laboratory.



The unique structure of our high performance core is designed for fast dispersion and encapsulation of urine, which results in a dry surface even after multiple voids.



Odor control is a critical factor in ensuring discreteness for the user. ABENA's incontinence products are designed and engineered to ensure that minimal Odor arises from the product.



Effective leakage barriers, both visible and invisible and latex-free elastics ensure a snug fit, which reduces the risk of leakages.



Wetness indicator is a liquid-sensitive strip that changes colour or fades out indicating when it is time to change the product.



A combination of 100% breathability, Top-Dry and leakage protection provides a dry surface and healthy skin even after multiple voids maintaining the skin's natural pH balance.



The product is fully breathable for the comfort of the user; also during long-time use. Helps maintain skin integrity.



Unique Top-Dry Technology that provides a rapid absorption and ensures a dry surface even after multiple voids.



The product does not contain latex.



Drop system indicating absorption. The more drops, the higher the absorption. Ranges from 1 to 9 drops.



The bag is 30% PCR, which means made from plastic materials used by consumers. This will allow us to reduce our plastic waste and to help take care of the environment.



ABENA's unique feature preventing leakage all around the product and providing an optimal fit.

