



Product Main Number 1999915563



Underpad, ABENA Washable Underpads, 90x85cm, with flaps

- ✓ High absorbency
- ✓ Washable material is temperature-safe up to 95 °C
- ✓ Peace of mind around the clock
- ✓ Soft and comfortable material



Product Description

ABENA Washable Underpads with flaps are primarily used as protective covers for beds and mattresses, providing user peace of mind in any situation. The underpad is particularly suitable for incontinent users who want added but discreet protection. The underpad is washable at temperatures of up to 95 °C. The underpad has an absorbent surface and a liquid safe back sheet, which, combined with the sealed edges, allow the pad to provide optimum leakage protection. In addition to helping provide peace of mind when sitting or lying down, the use of underpads also helps reduce the risk of contamination. At the same time, the underpad ensures increased hygiene. The surface is soft and comfortable, and the underpad is breathable, which ensures optimal comfort for the user and reduces the risk of skin irritation. At the same time, the product has rounded corners to reduce wear and tear from laundering. ABENA Washable with side flaps is specially designed for active and restless users. The side flaps ensure that the underpad does not move.



Specifications

| | |
|----------------------------------|--|
| Base name | Underpad |
| Brand | ABENA |
| Sub-Brand | Washable Underpads |
| Secondary color | Light blue |
| Properties | With flaps |
| Single or multiple use | Multiple use |
| Material | Cotton, polyester, rayon, TPU |
| Length/depth | 90 cm |
| Width | 85 cm |
| Absorption Capacity | 2500 ml |
| Certificates | CE. OEKO-TEX. MD. |
| CE Class (Medical Devices) | Class I |
| Product or test standards | ISO 11948-1 |
| Directives, regulations and acts | MDR (EU) 2017/745 |
| Shelf Life | 5 years |
| Storage Instructions | Store dry, room temperature and no direct sunlight. |
| Washing Instructions | Washable at 95 °C. Withstands tumble drying at medium temperature. Avoid fabric softener, chlorine and solvents, as this impairs the product's absorbency. Withstands at least 100 wash cycles. Product should be washed before use. |
| Product Disposal Instructions | Can be disposed of with normal household waste sorted according to local regulations. |
| Packaging Disposal Instructions | Can be recycled or incinerated. |



Instructions for use/application

The underpad must be placed on the mattress with the blue side down. Make sure the underpad is completely smoothed out. For optimum protection, the product should be washed before use.

Packaging data

| Unit | Contains | Length | Width | Height | EAN |
|------|----------|--------|---------|---------|---------------|
| cil | 30 pcs | 48 cm | 31.5 cm | 63.5 cm | 5713571007143 |



Regulation (EU) 2017/745 is a regulation of the European Union on the clinical investigation and sale of medical devices for human use. It repeals Directive 93/42/EEC, which concerns medical devices, and Directive 90/385/EEC, which concerns active implantable medical devices, on 26 May 2021.



The CE mark guarantees that a product is safe to use and complies with all safety precautions. CE stands for Conformité Européenne (European Conformity) and is mainly found on electronic equipment, safety equipment, construction products, and medical devices.



OEKO-TEX, also known as “Confidence in Textiles”, is a health and environmental label used worldwide. The certificate ensures that textiles do not contain substances that are harmful to health or the environment.



The product is fully breathable for the comfort of the user; also during long-time use. Helps maintain skin integrity.

